


Biotech SmartSaver

- Modern, compact design
- Powered from PC USB port
- LED indication of WASTE/RECYCLE positions
- Manual WASTE/RECYCLE control
- Compatible with any HPLC detector
- Easy Plug & Play Installation
- Drivers and Software Included

SmartSaver recycler saves up to 90% of mobile phase by redirection of the pure solvent to the solvent reservoir during the isocratic HPLC. SmartSaver brings quite new approach to solvent recycler design.

Very compact instrument is powered directly from chromatography data system PC, no power adapter is required. User friendly software is provided to configure the parameters, and to perform on-line monitoring/audit trail. Analog input allows unipolar or bipolar operation within range of $\pm 1V$. TTL/contact closure can be configured as start, auto-zero or valve position control input.

SmartSaver continuously monitors output signal of the chromatographic detector. Analog-to-digital converter produces digital data for further evaluation in build-in processor. If the input signal level exceeds the preset value (threshold), SmartSaver redirects the flow to the waste, respecting the transport time from detector to the switching valve. When signal goes down SmartSaver is waiting for transport delay again and then switch the mobile phase back to the reservoir. Autosampler injection marker connected to the SmartSaver may ensure the zeroing of signal input at the moment of injection



Specifications

- Powered from PC USB port
- Input range $\pm 1V$, optimized for "integrator" detector output
- 12 bit analog-to-digital converter
- 1 Hz data rate
- LED indication of WASTE/RECYCLE positions
- Manual WASTE/RECYCLE control
- User configurable TTL/contact closure input (start, valve position, zero ...)
- Compatible with any HPLC detector
- Wetted material: PEEK
- Connection: $\frac{1}{4}$ -28 Flat Bottom
- Maximum pressure: 30 psi/0.2 MPa
- Only two operating parameters – Threshold, Transport Delay
- Easy Plug&Play Installation
- Drivers and Software Included

SmartSaver comes with Standard Universal Cables, Instrument signal connections, BNC Coaxial connector for Agilent Detectors Output, LC connection tubing (PTFE & PEEK), LC connectors, LC Ferrules and Operations Manual with System CD

Ordering Information

Description	QTY	Part No.
Biotech SmartSaver		
Model SmartSaver	1	05AF-4011