


# Gas Sampling Bags

## Gas Sampling Bags with Swagelog® SS Fitting for Septum/Tubing



### Ordering Information

Description	QTY	Part No.
<b>ALTEF® Bags</b>		
6 x 6 inch, 0.6 L	1	061R-7101
7 x 7 inch, 1.0 L	1	061R-7102
6 x 10 inch, 1.2 L	1	061R-7103
6 x 12 inch, 1.5 L	1	061R-7104
9 x 9 inch, 2.0 L	1	061R-7105
12 x 12 inch, 5.0 L	1	061R-7106
12 x 17 inch, 8.0 L	1	061R-7107
12 x 19 inch, 10.0 L	1	061R-7108
18 x 18 inch, 16.0 L	1	061R-7109
18 x 24 inch, 25.0 L	1	061R-7110
24 x 24 inch, 40.0 L	1	061R-7111
24 x 30 inch, 56.0 L	1	061R-7112
24 x 36 inch, 73.0 L	1	061R-7113
30 x 30 inch, 80.0 L	1	061R-7114
30 x 36 inch, 100.0 L	1	061R-7115
<b>TEFLON® FEP Bags</b>		
6 x 6 inch, 0.6 L	1	061R-7201
7 x 7 inch, 1.0 L	1	061R-7202
6 x 10 inch, 1.2 L	1	061R-7203
6 x 12 inch, 1.5 L	1	061R-7204
9 x 9 inch, 2.0 L	1	061R-7205
10 x 10 inch, 3.0 L	1	061R-7206
12 x 12 inch, 5.0 L	1	061R-7207
12 x 18 inch, 9.5 L	1	061R-7208
12 x 19 inch, 10.0 L	1	061R-7209
18 x 18 inch, 16.0 L	1	061R-7210
24 x 24 inch, 40.0 L	1	061R-7211
24 x 30 inch, 56.0 L	1	061R-7212
24 x 36 inch, 73.0 L	1	061R-7213
30 x 30 inch, 80.0 L	1	061R-7214