

GC Manual Syringes

Micro-Volume Syringes

- Sample is contained in the needle with full plunger displacement, no dead volume
- High precision and accuracy down to 0.1µL
- Seal between plunger and needle can be tightened for high-pressure injection
- Replaceable needle/plunger kits and barrels



Ordering Information

Volume, Termination, Gauge, Length, PS#	QTY	OEM Model	OEM No.	Part No.
Hamilton Syringes				
0.5µL, KH, 25, 70mm, 2	1	7000.5	86259	24-86259
0.5µL, KH, 25, 70mm, 3	1	7000.5	86250	24-86250
1.0µL, KH, 22s, 70mm, 2	1	7101	86211	24-86211
1.0µL, KH, 25s, 70mm, 2	1	7001	80135	24-80135
1.0µL, KH, 25s, 70mm, 3	1	7001	80100	24-80100
2.0µL, KH, 25, 70mm, 2	1	7002	88411	24-88411
5.0µL, KH, 24, 70mm, 2	1	7105	88011	24-88011
5.0µL, KH, 24, 70mm, 3	1	7105	88000	24-88000
SGE Syringes				
0.5µL, R, 23, 70mm, 2	1	0.5BR-7BV	000311	21-000311
0.5µL, R, 23, 70mm, Cone	1	0.5BR-7	000310	21-000310
0.5µL, R, 26, 70mm, Cone	1	0.5BR-OC-7/0.47	000376	21-000376
1.0µL, R, 23, 70mm, 2	1	1BR-7BV	000506	21-000506
1.0µL, R, 23, 70mm, Cone	1	1BR-7	000505	21-000505
1.0µL, R, 23, 50mm, Cone	1	1BR-5	000500	21-000500
1.0µL, R, 26, 70mm, Cone	1	1BR-7/0.47	000570	21-000570
5.0µL, R, 23, 70mm, 2	1	5BR-7BV	000803	21-000803
5.0µL, R, 23, 70mm, Cone	1	5BR-7	000802	21-000802
5.0µL, R, 23, 50mm, Cone	1	5BR-5	000800	21-000800
SGE Replacement Plunger-in-Needle Kits				
0.5µL, 23, 70mm, 2	1	NP0.5B-7BV	033060	21-033060
0.5µL, 23, 70mm, Cone	1	NP0.5B-7C	033057	21-033057
0.5µL, 26, 70mm, Cone	1	NP0.5B-OC-7	033630	21-033630
1.0µL, 23, 70mm, 2	1	NP1B-7BV	034060	21-034060
1.0µL, 23, 70mm, Cone	1	NP1B-7C	034057	21-034057
1.0µL, 23, 50mm, Cone	1	NP1B-5C	034055	21-034055
1.0µL, 26, 70mm, Cone	1	NP1B-OC-7/0.47	034610	21-034610
5.0µL, 23, 70mm, 2	1	NP5B-7BV	035058	21-035058
5.0µL, 23, 70mm, Cone	1	NP5B-7C	035057	21-035057
5.0µL, 23, 50mm, Cone	1	NP5B-5C	035055	21-035055