

## Gas-Tight Syringes with Shut-Off Valves

- Ideal for sample storage and transportation
- Push-pull shut-off valve
- Use VICI A-2 for VOC gas sampling from a high-pressure cylinder through a high purity VOC single stage regulator



| Ordering Information                               |     |                |         |           |
|----------------------------------------------------|-----|----------------|---------|-----------|
| Syringe Volume                                     | QTY | OEM Model      | OEM No. | Part No.  |
| <b>VICI A-2 Syringe with Push Botton Valve</b>     |     |                |         |           |
| 25µL                                               | 1   |                | 050023  | 25-050023 |
| 50µL                                               | 1   |                | 050024  | 25-050024 |
| 100µL                                              | 1   |                | 050025  | 25-050025 |
| 250µL                                              | 1   |                | 050031  | 25-050031 |
| 500µL                                              | 1   |                | 050032  | 25-050032 |
| 1mL                                                | 1   |                | 050033  | 25-050033 |
| 2mL                                                | 1   |                | 050034  | 25-050034 |
| 5mL                                                | 1   |                | 050035  | 25-050035 |
| 10mL                                               | 1   |                | 050036  | 25-050036 |
| <b>SGE Syringe with Luer-Lock Valve*</b>           |     |                |         |           |
| 1mL                                                | 1   | 1MR-VLL-GT     | 008160  | 21-008160 |
| 2.5mL                                              | 1   | 2.5MDR-VLL-GT  | 008560  | 21-008560 |
| 5mL                                                | 1   | 5MDR-VLLMA-GT  | 008770  | 21-008770 |
| 10mL                                               | 1   | 10MDR-VLLMA-GT | 008970  | 21-008970 |
| 25mL                                               | 1   | 25MDR-VLLMA-GT | 009472  | 21-009472 |
| 50mL                                               | 1   | 50MR-VLLMA-GT  | 009670  | 21-009670 |
| 100mL                                              | 1   | 100MR-VLLMA-GT | 009770  | 21-009770 |
| <b>SGE Syringe with Removable Needle and Valve</b> |     |                |         |           |
| 50µL                                               | 1   | 50R-V-GT       | 004279  | 21-004279 |
| 100µL                                              | 1   | 100R-V-GT      | 005279  | 21-005279 |
| 250µL                                              | 1   | 250R-V-GT      | 006279  | 21-006279 |
| 500µL                                              | 1   | 500R-V-GT      | 007279  | 21-007279 |
| 1mL                                                | 1   | 1MDR-V-GT      | 008110  | 21-008110 |
| 2.5mL                                              | 1   | 2.5MDR-V-GT    | 008510  | 21-008510 |

\* Needles sold separately. See next page.

## Needles for Syringes with Shut-Off Valves

| Ordering Information                                           |     |           |         |           |
|----------------------------------------------------------------|-----|-----------|---------|-----------|
| Syringe Volume, Gauge, Length, PS#                             | QTY | OEM Model | OEM No. | Part No.  |
| <b>Needles for VICI A-2 Syringe with Push Botton Valve</b>     |     |           |         |           |
| 1mL–100mL, 26, 50mm, 2                                         | 3   |           | 943050  | 25-943050 |
| 1mL–100mL, 26s, 50mm, 2                                        | 3   |           | 943051  | 25-943051 |
| <b>Needles for SGE Syringe with Luer-Lock Valve</b>            |     |           |         |           |
| 1mL–100mL, 23, 50mm, 2                                         | 5   | NLL-5/23  | 039802  | 21-039802 |
| 1mL–100mL, 18, 50mm, 2                                         | 5   | NLL-5/18  | 039842  | 21-039842 |
| 1mL–100mL, 22, 51mm, 3                                         | 2   | NLL-LC    | 039895  | 21-039895 |
| <b>Needles for SGE Syringe with Removable Needle and Valve</b> |     |           |         |           |
| 25µL–2.5mL, 23, 50mm, Cone                                     | 1   |           | 031905  | 21-031905 |