



Vespel ferrules

Ordering Information			
Description	QTY	Agilent No.	Part No.
Vespel Ferrules			
1/16" to 0.4mm, short, for 0.25mmID columns	10	5181-3322	062M-4016
1/16" to 0.5mm, short, for 0.32mmID columns	10	5062-3513	062M-4026
1/16" to 0.8mm, short, for 0.53mmID columns	10	5062-3511	062M-4036
1/16" to 0.4mm, for 0.25mmID columns	10	5182-0740	062M-4046
1/16" to 0.5mm, for 0.32mmID columns	10	5182-0741	062M-4056
1/16" to 0.8mm, for 0.53mmID columns	10	5182-0742	062M-4066
1/16" 2-hold, 0.5mm, for 0.32mmID columns	10		062M-4076
1/16" ID Straight	10		062M-4086
1/8" ID Straight	10		062M-4096
1/4" ID Straight	10	0700-1327	062M-4106
1/8" Reducing to 1/16" for 1/16"OD columns	10		062M-4116
1/4" Reducing to 1/16", for 1/16"OD tubing	10		062M-4126
1/4" Reducing to 1/8", for 1/8"OD tubing	10		062M-4136
1/16" Blank, no holes	10	5021-7133	062M-4146
1/8" Blank, no holes	10		062M-4156
1/4" Blank, no holes	10		062M-4166

Capillary Nuts



Capillary column nuts

Ordering Information			
Description	QTY	Agilent No.	Part No.
Capillary Nuts			
Stainless steel capillary nut for short ferrules	2	5181-8830	062M-5012
Stainless steel capillary nut for standard ferrules	2	05921-21170	062M-5022

SilTite Ferrules

- Metal ferrules for fused silica GC columns and tubing
- Ideal for GC/MS
- Eliminate contamination from vespel or graphite
- Eliminate annoying leaks
- Reduce air background
- No more retightening of Vespel ferrules after temperature cycling
- Reduce GC and MS down time



SilTite nuts and ferrules

Ordering Information			
Description	QTY	Agilent No.	Part No.
SilTite Ferrules Kit : 10 Ferrules, 2 Nuts			
For MSD, 0.10 to 0.25mmID Columns	1	5181-8361	062M-6011
For MSD, 0.32mmID Columns	1	5181-8362	062M-6021
For MSD, 0.45 to 0.53mmID Columns	1	5181-8363	062M-6031
For Split/Spliless, 0.10 to 0.25mmID Columns	1		062M-6041
For Split/Spliless, 0.32mmID Columns	1		062M-6051
For Split/Spliless, 0.45 to 0.53mmID Columns	1		062M-6061
SilTite Ferrules Kit : 10 Ferrules, 2 Nuts, 2 Inlet Base Seals			
For Split/Spliless, 0.10 to 0.25mmID Columns	1		062M-6071
For Split/Spliless, 0.32mmID Columns	1		062M-6081
For Split/Spliless, 0.45 to 0.53mmID Columns	1		062M-6091
Replacement			
SilTite Ferrules 0.4mm for 0.1 to 0.25mmID Columns	10	5188-5361	062M-6106
SilTite Ferrules 0.5mm for 0.32mmID Columns	10	5188-5362	062M-6116
SilTite Ferrules 0.8mm for 0.53mmID Columns	10	5188-5363	062M-6126
SilTite Nuts	5		062M-6135