

PerkinElmer Supplies/Parts

Capillary Adaptors



Ordering Information			
Description	QTY	PerkinElmer No.	Part No.
SilTite Ferrules Kit : 10 Ferrules, 2 Nuts			
For detectors, PerkinElmer style	1	N6120020	064M-8011
For detectors, Swagelog style	1		064M-8021
For injectors, PerkinElmer style	1	N6100157	064M-8031
For injectors, Swagelog style	1		064M-8041

SilTite Ferrules

- Metal ferrules for fused silica GC columns and tubing
- Ideal for GC/MS
- Eliminate contamination from vespel or graphite
- Eliminate annoying leaks
- Reduce air background
- No more retightening of Vespel ferrules after temperature cycling
- Reduce GC and MS down time



SilTite nuts and ferrules

Ordering Information			
Description	QTY	PerkinElmer No.	Part No.
SilTite Ferrules Kit : 10 Ferrules, 2 Nuts			
For MSD, 0.10 to 0.25mmID Columns	1		064M-6011
For MSD, 0.32mmID Columns	1		064M-6021
For MSD, 0.45 to 0.53mmID Columns	1		064M-6031
Replacement			
SilTite Ferrules 0.4mm for 0.1 to 0.25mmID Columns	10		064M-6046
SilTite Ferrules 0.5mm for 0.32mmID Columns	10		064M-6056
SilTite Ferrules 0.8mm for 0.53mmID Columns	10		064M-6066
SilTite Nuts	5		064M-6075