

Varian Supplies/Parts

Autosampler Syringes



063A-1011 063A-2031




063A-1056

Needle Termination Codes	
Hamilton	
N	Cemented Needle
RN	Removable Needle
SGE	
F	Fixed Needle
R	Removable Needle
GT	Gas-Tight

Ordering Information

Volume, Code, Gauge, Length, Point Style	QTY	OEM No.	Varian No.	Part No.
Hamilton Syringes				
10µL, RN, 26s, 50mm, #5 for 8100/8200	1	202880	391898600	063A-1011
5µL, N, 26s, 50mm, #2 for 8400/9000	1	87900	HM87900	063A-1021
10µL, N, 26s, 50mm, #2 for 8400/9000	1	80300	HM80300	063A-1031
10µL, N, 26s, 50mm, #2 for 8400/9000	6	80366	HM80366	063A-1046
Needles 5-10µL, 23s, 50mm, #5	6	202903	391898700	063A-1056
SGE Syringes				
10µL, F, 26s, GT, 50mm, #5 for 8100/8200	1	002923		063A-2011
10µL, R, GT, 26s, 50mm, #5 for 8100/8200	1	002924		063A-2021
100µL, R, 26s, 50mm, #5 for CP 8100/8200	1	005921	391898600	063A-2031
10µL, F, 26s, 50mm, Cone, for CP 8400/9000	1	001800		063A-2041
10µL, F, 26s, 50mm, #2, for CP 8400/9000	1	001804	HM80300	063A-2051
Needles 10µL, 23s, 50mm, #5	2	037779	391898700	063A-2062
Needles 10µL, 26s, 50mm, #5	1	037780	391898700	063A-2071
Needles 10µL, 26s, 50mm, #2	2	037776		063A-2082



Point Style #5
Side-hole cone tip



Point Style #2
Beveled needle tip



Point Style Cone
Cone tip