

Varian Supplies/Parts



Teflon ferrules

Ordering Information			
Description	QTY	Varian No.	Part No.
Teflon Ferrules			
1/16" to 0.3mm, for 0.18mm and smaller ID columns	10		063M-7016
1/16" to 0.425mm, for 0.25mmID columns	10		063M-7026
1/16" to 0.5mm, for 0.32mmID columns	10		063M-7036
1/16" to 0.8mm, for 0.53mmID columns	10		063M-7046
1/16" to 1.0mm, for 0.65mmID columns	10		063M-7056
1/16" to 1.2mm, for 0.75mmID columns	10		063M-7066
1/16" Straight	10		063M-7076
1/8" Straight	10		063M-7086
3/16" Straight	10		063M-7096
1/4" Straight	10	CR214400	063M-7106
1/8" Reducing to 0.425mm, for 0.25mmID columns	10		063M-7116
1/8" Reducing to 0.5mm, for 0.32mmID columns	10		063M-7126
1/8" Reducing to 0.8mm, for 0.53mmID columns	10		063M-7136
1/8" Reducing to 1/16" for 1/16"OD columns	10		063M-7146
1/8" Reducing to 2-hole, 0.5mm, for 0.32mmID columns	10		063M-7156
1/4" Reducing to 6mm, for 1/16"OD tubing	10		063M-7166
1/4" Reducing to 1/16", for 1/16"OD tubing	10		063M-7176
1/4" Reducing to 1/8", for 1/8"OD tubing	10		063M-7186
1/4" Reducing to 3/16", for 3/16"OD tubing	10		063M-7196
1/16" Blank, no holes	10		063M-7206
1/8" Blank, no holes	10		063M-7216
1/4" Blank, no holes	10		063M-7226

Capillary Nuts



Capillary Nuts

Ordering Information			
Description	QTY	Varian No.	Part No.
Capillary Nuts			
Brass Capillary Nut	2	394955100	063M-5012
Stainless Steel Capillary Nut	2	CP743117	063M-5022

SilTite Ferrules

- Metal ferrules for fused silica GC columns and tubing
- Ideal for GC/MS
- Eliminate contamination from vespel or graphite
- Eliminate annoying leaks
- Reduce air background
- No more retightening of Vespel ferrules after temperature cycling
- Reduce GC and MS down time



SilTite nuts and ferrules

Ordering Information			
Description	QTY	Varian No.	Part No.
SilTite Ferrules Kit : 10 Ferrules, 2 Nuts			
For 0.10 to 0.25mmID Columns	1	SG073300	063M-6011
For 0.32mmID Columns	1	SG073301	063M-6021
For 0.45 to 0.53mmID Columns	1	SG073302	063M-6031
Replacement			
SilTite Ferrules 0.4mm for 0.1 to 0.25mmID Columns	10	SG073220	063M-6046
SilTite Ferrules 0.5mm for 0.32mmID Columns	10	SG073221	063M-6056
SilTite Nuts	5	SG073224	063M-6065