



Short ferrules



Straight ferrules



Reducing ferrules



VespeL ferrules

Ordering Information			
Description	QTY	Varian No.	Part No.
40% Graphite/ 60% VespeL Ferrules			
1/16" to 0.3mm, for 0.18mm and smaller ID columns	10	CR213103	063M-3016
1/16" to 0.425mm, for 0.25mmID columns	10	CR213104	063M-3026
1/16" to 0.5mm, for 0.32mmID columns	10	CR213105	063M-3036
1/16" to 0.8mm, for 0.53mmID columns	10	CR213108	063M-3046
1/16" to 1.0mm, for 0.65mmID columns	10	CR213110	063M-3056
1/16" to 1.2mm, for 0.75mmID columns	10	CR213112	063M-3066
1/16" 2-hole, 0.425mm, for 0.25mmID columns	10	CR213124	063M-3076
1/16" 2-hole, 0.5mm, for 0.32mmID columns	10	CR213125	063M-3086
1/16" to 0.425mm Short, for 0.25mmID columns	10	CR213164	063M-3096
1/16" to 0.5mm Short, for 0.32mmID columns	10	CR213165	063M-3106
1/16" to 0.8mm Short, for 0.53mmID columns	10	CR213168	063M-3116
1/16" Straight	10		063M-3126
1/8" Straight	10	CR213200	063M-3136
1/4" Straight	10	CR213400	063M-3146
1/8" Reducing to 0.425mm, for 0.25mmID columns	10	CR213204	063M-3156
1/8" Reducing to 0.5mm, for 0.32mmID columns	10	CR213205	063M-3166
1/8" Reducing to 0.8mm, for 0.53mmID columns	10	CR213208	063M-3176
1/8" Reducing to 1/16" for 1/16"OD columns	10		063M-3186
1/8" Reducing to 2-hole, 0.5mm, for 0.32mmID columns	10		063M-3196
1/4" Reducing to 0.5mm, for 0.32mmID columns	10	CR213405	063M-3206
1/4" Reducing to 6mm, for 1/16"OD tubing	10		063M-3216
1/4" Reducing to 0.8mm, for 0.53mmID columns	10		063M-3226
1/4" Reducing to 1/16", for 1/16"OD tubing	10		063M-3236
1/4" Reducing to 1/8", for 1/8"OD tubing	10		063M-3246
1/16" Blank, no holes	10	CR2131XX	063M-3256
1/8" Blank, no holes	10		063M-3266
1/4" Blank, no holes	10		063M-3276
VespeL Ferrules			
1/16" to 0.3mm, for 0.18mm and smaller ID columns	10	CR212103	063M-4016
1/16" to 0.425mm, for 0.25mmID columns	10	CR212104	063M-4026
1/16" to 0.5mm, for 0.32mmID columns	10	CR212105	063M-4036
1/16" to 0.8mm, for 0.53mmID columns	10	CR212108	063M-4046
1/16" to 1.0mm, for 0.65mmID columns	10		063M-4056
1/16" to 1.2mm, for 0.75mmID columns	10		063M-4066
1/16" 2-hole, 0.425mm, for 0.25mmID columns	10		063M-4076
1/16" 2-hole, 0.5mm, for 0.32mmID columns	10	CR212125	063M-4086
1/16" Straight	10		063M-4096
1/8" Straight	10		063M-4106
1/4" Straight	10		063M-4116
1/8" Reducing to 0.425mm, for 0.25mmID columns	10		063M-4126
1/8" Reducing to 0.5mm, for 0.32mmID columns	10		063M-4136
1/8" Reducing to 0.8mm, for 0.53mmID columns	10		063M-4146
1/8" Reducing to 1/16" for 1/16"OD columns	10		063M-4156
1/8" Reducing to 2-hole, 0.5mm, for 0.32mmID columns	10		063M-4166
1/4" Reducing to 0.5mm, for 0.32mmID columns	10		063M-4176
1/4" Reducing to 6mm, for 1/16"OD tubing	10		063M-4186
1/4" Reducing to 0.8mm, for 0.53mmID columns	10		063M-4196
1/4" Reducing to 1/16", for 1/16"OD tubing	10		063M-4206
1/4" Reducing to 1/8", for 1/8"OD tubing	10		063M-4216
1/16" Blank, no holes	10		063M-4226
1/8" Blank, no holes	10		063M-4236
1/4" Blank, no holes	10		063M-4246