

SPE Bulk Packings and Accessories

SPE Empty Tube and Frits



VertiPak™ empty tube are polypropylene with male luer outlets. Frits are 20µm polyethylene. Select empty tube and loose frits to pack your own custom SPE tubes. Empty tubes with 2 bottom frits installed work well as filter tubes to remove particulates down to 20µm.

Empty Tube and Frits		
Description	QTY	Part No.
Empty Tube		
1 mL	100	0133-0001
1 mL	500	0133-0002
3 mL	100	0133-0003
6 mL	100	0133-0004
10 mL	100	0133-0005
12 mL	100	0133-0006
35 mL	100	0133-0007
60 mL	50	0133-0008
20µm Loose Frits		
For 1 mL Reservoir	100	0133-0009
For 3 and 10 mL Reservoir	100	0133-0010
For 6 mL Reservoir	100	0133-0011
For 12 mL Reservoir	100	0133-0012
For 35 mL Reservoir	100	0133-0013
For 60 mL Reservoir	100	0133-0014

Custom Packed SPE Packings

Vertical will custom request on your letterhead for a price quotation. Be sure to include your contact.

Please specify:

1. Packing materials
2. Quantity
3. Format size and style

Custom orders are not returnable. A minimum order may be required.

SPE Caps and Adaptors



SPE Caps and Adaptors

Description	QTY	Part No.
Inlet Caps for Cartridges	50	0133-0019
Inlet Caps for 1 mL	50	0133-0020
Inlet Caps for 3 mL	50	0133-0021
Inlet Caps for 6 mL	50	0133-0022
Inlet Caps for 12 mL	50	0133-0023
Inlet Caps for 35 mL	50	0133-0024
Inlet Caps for 60 mL	50	0133-0025
Outlet Caps for Male Luers	50	0133-0026
Adaptors for 1,3,6,12mL Tubes	5	0133-0027
Adaptors for 35,60 mL Tubes	5	0133-0028

Bulk SPE Packings

Bulk SPE Packings

Description	QTY	Part No.
C18 (17%)	100g	0133-0031
C8	100g	0133-0032
C2	100g	0133-0033
Phenyl	100g	0133-0034
Silica	100g	0133-0035
NH2	100g	0133-0036
PSA	100g	0133-0050
CN	100g	0133-0044
Diol	100g	0133-0037
Florisil	100g	0133-0038
Florisil-PR	100g	0133-0045
AL-N	100g	0133-0046
AL-A	100g	0133-0047
AL-B	100g	0133-0048
SCX	100g	0133-0039
SCX-2	100g	0133-0051
CBA	100g	0133-0049
SAX	100g	0133-0040
Carbograph	100g	0133-0041
C8/SCX	100g	0133-0042
C8/SAX	100g	0133-0043