

Vacuum Manifolds and Accessories



Vacuum Manifold Components

- A. Glass Chamber
- B. Vacuum Valve and Gauge
- C. Polypropylene Lid
- D. Stopcock Valves
- E. Collection Rack Plates
- F. Support Posts for Collection Racks
- G. Retaining Clips for Collection Racks
- H. Lid Legs
- I. Manifold Inlet Caps
- J. Polypropylene Needles
- K. SPE Reservoirs (sold separately)

- 12- and 24-Port Manifolds
- Glass Chamber for Visual Monitoring
- Accepts Standard Male Luer Devices

Vacuum manifolds allow you to process multiple samples simultaneously, saving time and effort. Manifold systems come complete with the components listed below. Stainless steel or Teflon® needles are available separately.

Includes Glass Chamber, Vacuum Valve and Gauge, Polypropylene Lid, Stopcock Valves, Collection Racks, Support Posts, Retaining Clips, Manifold Inlet Caps, Lid Legs, and Polypropylene Needles.

The 12-Port Manifold also includes a waste container.



12-Port Manifold



24-Port Manifold

12-Port Manifold

- For up to 12 Samples
- Includes one waste container

12-Port Vacuum Manifold		
Description	QTY	Part No.
12-Port Vacuum Manifold	Ea	0134-0100
Replacement Parts		
Lid, Gaskets, and 12 Stopcocks	Ea	0134-0102
Glass Chamber	Ea	0134-0103
Vacuum Gauge, Valve, and Glass Chamber	Ea	0134-0104
Collection Rack, 12-Port Size*	Ea	0134-0105
Gaskets, 12-Port Size	2	0134-0106
One-Way Stopcocks	12	0134-0107
Waste Container	2	0134-0108

*12-Port Collection Rack consists of 3 support posts, bottom plate, 13- and 16mm plates, autosampler vial plate, volumetric plate, and 12 retaining clips

24-Port Manifold

- Process up to 24 Samples

24-Port Manifold & Accessories		
Description	QTY	Part No.
24-Port Vacuum Manifold	Ea	0134-0300
Replacement Parts		
Lid, Gaskets, and 24 Stopcocks	Ea	0134-0301
Glass Chamber	Ea	0134-0302
Vacuum Gauge, Valve, and Glass Chamber	Ea	0134-0303
Collection Rack, 24-Port Size**	Ea	0134-0304
Gaskets, 24-Port Size	2	0134-0205
One-Way Stopcocks	24	0134-0306

**24-Port Collection Racks consist of 3 support posts, bottom plate, dimple plate, 13- and 16mm plates, and 12 retaining clips.

SPE