

Membrane Filters and Accessories

Membrane Filters

- Variety of Pore Sizes and Membrane Materials
- Non-Sterile



Hydrophilic : Aqueous and organic solvent

Ordering Information			
Size	Pore Size	QTY	Part No.
NYLON Membrane Filters			
13mm	0.2µm	100	0231-0101
	0.45µm	100	0231-0301
25mm	0.2µm	100	0233-0101
	0.45µm	100	0233-0301
47mm	0.2µm	100	0235-0101
	0.45µm	100	0235-0301
CA Membrane Filters			
13mm	0.2µm	100	0231-1101
	0.45µm	100	0231-1301
25mm	0.2µm	100	0233-1101
	0.45µm	100	0233-1301
47mm	0.2µm	100	0235-1101
	0.45µm	100	0235-1301
MCE Membrane Filters			
13mm	0.2µm	100	0231-8101
	0.45µm	100	0231-8301
25mm	0.2µm	100	0233-8101
	0.45 µm	100	0233-8301
47mm	0.2 µm	100	0235-8101
	0.45 µm	100	0235-8301
PVDF-HL Membrane Filters			
13mm	0.2µm	100	0231-3101
	0.45µm	100	0231-3301
25mm	0.2µm	100	0233-3101
	0.45µm	100	0233-3301
47mm	0.2µm	100	0235-3101
	0.45µm	100	0235-3301

Hydrophobic : Non-Aqueous, organic and high solvent resistant

Ordering Information			
Size	Pore Size	QTY	Part No.
PTFE Membrane Filters			
13mm	0.2µm	100	0231-2101
	0.45µm	100	0231-2301
25mm	0.2µm	100	0233-2101
	0.45µm	100	0233-2301
47mm	0.2µm	100	0235-2101
	0.45µm	100	0235-2301
PVDF-HB Membrane Filters			
13mm	0.2µm	100	0231-9101
	0.45µm	100	0231-9301
25mm	0.2µm	100	0233-9101
	0.45µm	100	0233-9301
47mm	0.2µm	100	0235-9101
	0.45µm	100	0235-9301

Syringe Filter Holder

Use this reusable syringe filter housing with disposable filter membranes to filter samples or small quantities of mobile phases. The stainless steel or polypropylene filter housing holds your choice of membrane. The outlet accepts any standard luer-hub needle.



Ordering Information		
Description	QTY	Part No.
Syringe Filter Holder		
13mm Stainless Steel	1	0210-0001
25mm Stainless Steel	1	0210-0002
13mm Polypropylene Steel	10	0210-0003
25mm Polypropylene Steel	12	0210-0004


Membrane Filters and Accessories

Glass Syringes

- Reinforced flanges with two flat sides to prevent rolling
- All borosilicate glass plunger and barrel for accurate, leak-proof performance and resistance to chemicals and thermal shock.
- Graduations are fused on to the syringe barrel
- Interchangeable plungers



5mL Glass syringe with glass luer slip tip



20mL Glass syringe with metal luer lock tip

Plastic Syringes

- Made of medical-grade polypropylene
- Sterile Disposable
- Ideal for use with syringe filters and general purpose applications
- Come in an easy-to-open sealed package
- Eliminate Ion Contamination in IC Applications



Ordering Information		
Size	QTY	Part No.
Glass Syringe, Luer-Tip		
2mL	1	0211-1211
	12	0211-1213
5mL	1	0211-1411
	12	0211-1413
10mL	1	0211-1511
	12	0211-1513
20mL	1	0211-1611
	12	0211-1613
50mL	1	0211-1811
	12	0211-1813
Glass Syringe, Luer-Lock		
2mL	1	0211-2211
	12	0211-2213
5mL	1	0211-2411
	12	0211-2413
10mL	1	0211-2511
	12	0211-2513
20mL	1	0211-2611
	12	0211-2613
50mL	1	0211-2811
	12	0211-2813

Ordering Information		
Size	QTY	Part No.
Plastic Syringe, Luer-Tip		
1mL	100	0211-1125
3mL	100	0211-1325
5mL	100	0211-1425
10mL	100	0211-1525
20mL	50	0211-1624
Plastic Syringe, Luer-Lock		
3mL	100	0211-2325
5mL	100	0211-2425
10mL	100	0211-2525
20mL	50	0211-2624