


Spiral PEEK Tubings

- Self-adjusting length
- Will not bend or kink
- Spiral tubing expands and contracts and is made of PEEK polymer, a biocompatible, chemically inert material
- The coiled format springs back to keep excess tubing out of your way



Ordering Information		
Description	QTY	Part No.
Spiral Tubing, 1/16" OD x 0.005" ID, Red		
20cm, 0.47" Coil span	1	05AD-1361
40cm, 2.4" Coil span	1	05AD-1371
50cm, 3.0" Coil span	1	05AD-1381
100cm, 7.0" Coil span	1	05AD-1391
200cm, 13.0" Coil span	1	05AD-1401
Spiral Tubing, 1/16" OD x 0.010" ID, Blue		
20cm, 0.47" Coil span	1	05AD-1411
40cm, 2.4" Coil span	1	05AD-1421
50cm, 3.0" Coil span	1	05AD-1431
100cm, 7.0" Coil span	1	05AD-1441
200cm, 13.0" Coil span	1	05AD-1451
Spiral Tubing, 1/16" OD x 0.020" ID, Orange		
20cm, 0.47" Coil span	1	05AD-1461
40cm, 2.4" Coil span	1	05AD-1471
50cm, 3.0" Coil span	1	05AD-1481
100cm, 7.0" Coil span	1	05AD-1491
200cm, 13.0" Coil span	1	05AD-1501
Spiral Tubing, 1/16" OD x 0.030" ID, Green		
20cm, 0.47" Coil span	1	05AD-1511
40cm, 2.4" Coil span	1	05AD-1521
50cm, 3.0" Coil span	1	05AD-1531
100cm, 7.0" Coil span	1	05AD-1541
200cm, 13.0" Coil span	1	05AD-1551

SGE PEEKsil™ Tubings

- Ideal for LC systems where conventional 1/16" OD fittings are used
- Constructed from a core of smooth precision bore fused silica tubing bonded to a 1/16" OD PEEK™ outer sheath
- All PEEKsil tubing is pressure rated to 8,500 psi
- Resistant to strong acids and has an effective pH range of 0-10
- Not compatible with hydrofluoric acid



Ordering Information		
Description	QTY	Part No.
PEEKsil™ Tubing, 1/16" OD 25uL, Orange		
10cm	5	05AD-1565
15cm	5	05AD-1575
25cm	5	05AD-1585
50cm	2	05AD-1592
PEEKsil™ Tubing, 1/16" OD 50uL, Natural		
10cm	5	05AD-1605
15cm	5	05AD-1615
25cm	5	05AD-1625
50cm	2	05AD-1632
PEEKsil™ Tubing, 1/16" OD 75uL, Black		
10cm	5	05AD-1645
15cm	5	05AD-1655
25cm	5	05AD-1665
50cm	2	05AD-1672
PEEKsil™ Tubing, 1/16" OD 150uL, Red		
10cm	5	05AD-1685
15cm	5	05AD-1695
25cm	5	05AD-1705
50cm	2	05AD-1712
PEEKsil™ Tubing, 1/16" OD 175uL, Yellow		
10cm	5	05AD-1725
15cm	5	05AD-1735
25cm	5	05AD-1745
50cm	2	05AD-1782
PEEKsil™ Tubing, 1/16" OD 200uL, Blue		
10cm	5	05AD-1795
15cm	5	05AD-1775
25cm	5	05AD-1785
50cm	2	05AD-1792