

Valves Replacement Parts

Rheodyne® Rhebuild® Kit for Injection valves




RheBuild® Kits are a complete selection of all the parts you need for general valve maintenance. All are genuine Rheodyne® parts selected for your specific valve.

Kits for Standard Front-loading Injection Valves contain: Standard Rotor Seal, Stator Face Assembly, Isolation Seal, Needle Guide, Needle Port Cleaner, Hex Keys and Instructions. Kits for Models 7520 and 7526 contain: Inlet Stator, Rotor, Needle Guide, Needle Port Cleaner, Hex Key, Rotor Removal Tool and Instructions.

Kits for Rear-Loading Injection Valves contain: Standard Rotor Seal, Isolation Seal, Hex Keys, and Instructions.

Ordering Information		
Description	QTY	Part No.
Rheodyne® Rhebuilt Kit		
for 3725, 3725i, 3725-038, 3725i-038	1	05AE-2011
Kit to include stator for 7010	1	05AE-2021
for 7010, 7000	1	05AE-2031
for 7125, 7126	1	05AE-2041
for 7125-081	1	05AE-2051
for 7410	1	05AE-2061
for 7520, 7526	1	05AE-2071
for 7725, 7725i, 7726	1	05AE-2081
for 7750TPMV series	1	05AE-2091
for 8125, 8126	1	05AE-2101
for 9010, 9040	1	05AE-2111
for 9125, 9126	1	05AE-2121
for 9725, 9725i, 9726	1	05AE-2131
for 9750TPMV series	1	05AE-2141

Rheodyne® Replacement Parts



Ordering Information		
Description	QTY	Part No.
Rotor Seals		
Tefzel, for 7010, 7010-087, 7000, 7040	1	05AE-2151
Tefzel, for 7030	1	05AE-2161
Tefzel, for 7060, 7066, 9060	1	05AE-2171
Tefzel, for 7125, 7725	1	05AE-2181
Tefzel, for 7410	1	05AE-2191
Tefzel, for 7413	1	05AE-2201
Tefzel, for 8125	1	05AE-2211
Tefzel, for 9010	1	05AE-2221
Tefzel, for 9125, 9725	1	05AE-2231
PEEK, for 3725, 3725-038	1	05AE-2241
PEEK, for 7610-400, 7610-600	1	05AE-2251
PEEK, for 7125, 7725, 9125, 9725	1	05AE-2261
Stators		
for 3725, 3710-038, PEEK	1	05AE-2271
for 3725-038, 3710-038	1	05AE-2281
for 7125, 7010, 7000, 7030, 7040	1	05AE-2291
for 7125-081, 7010-087	1	05AE-2301
for 7060, 7066	1	05AE-2311
for 7410, 7413	1	05AE-2321
for 7610-600	1	05AE-2331
for 7725	1	05AE-2341
for 8125, 8126	1	05AE-2351
for 9125, 9010, 9030, 9725	1	05AE-2361
for 9060	1	05AE-2371
Stator Face Assemblies		
for 7125, 7126	1	05AE-2381
for 7725	1	05AE-2391
for 8125, 8126	1	05AE-2401
for 9125, 9010, 9030	1	05AE-2411
for 9060	1	05AE-2421
Miscellaneous Replacement Parts		
Needle guide (for front loading injectors)	1	05AE-2431
Isolation seal	1	05AE-2441