

HPLC Instrument Parts

Agilent HPLC 1100/1200

- Made of flexible stainless steel with 1/16" fittings for both ends
- Pre-swaged (p-s) fittings are assembled according to Swagelok specifications
- Non-swaged (n-s) capillaries include fittings, but are not assembled



05BH-0651

Ordering Information

Description	QTY	Agilent No.	Part No.
Universal Connecting Capillaries			
130mm x 0.25mmID, p-s to p-s, Pump to Autosampler	1	01090-87308	05BH-0631
160mm x 0.25mmID, p-s to p-s	1	G1313-87301	05BH-0641
90mm x 0.17mmID, n-s to n-s	1	G1316-87300	05BH-0651
130mm x 0.17mmID, p-s to n-s	1	01090-87305	05BH-0661
180mm x 0.17mmID, p-s to n-s Autosampler to Column	1	G1313-87305	05BH-0671
280mm x 0.17mmID, p-s to n-s Autosampler to Column	1	01090-87304	05BH-0681
380mm x 0.17mmID, p-s to n-s Column to DAD	1	G1315-87311	05BH-0691
500mm x 0.17mmID, n-s to n-s	1	G1328-87600	05BH-0701
600mm x 0.17mmID, p-s to n-s Pump to Autosampler	1	G1312-67305	05BH-0711
800mm x 0.17mmID, p-s to n-s Autosampler to Column	1	01048-87302	05BH-0721
35mm x 0.12mmID, p-s to p-s	1	79841-87609*	05BH-0731
70mm x 0.12mmID, n-s to n-s Column to Universal	1	G1316-87303	05BH-0741
105mm x 0.12mmID, p-s to n-s	1	01090-87611	05BH-0751
280mm x 0.12mmID, p-s to n-s Autosampler to Column	1	01090-87610	05BH-0761

*Made of 1.58mm (rigid) OD tubing

Beckman/Altex HPLC Systems



05CH-0011



05CH-0021



05CH-0081



05CH-0154



05CH-0181



05CH-0161



05CH-0171



05CH-0061



05CH-0041



05CH-0031



05CH-0071

Ordering Information

Description	QTY	Beckman No.	Part No.
Replacement parts for Beckman HPLC Systems			
Inlet Check Valve Assembly	1	240720	05CH-0011
Outlet Check Valve Assembly	1	240721	05CH-0021
Inlet Check Valve Cartridge for Check Valve Assembly	1	240620	05CH-0031
Outlet Check Valve Cartridge for Check Valve Assembly	1	240621	05CH-0041
Guide Sleeve for Graphite Plunger Guide	1	243713	05CH-0051
Outlet Check Filter Frit Assembly	1	240619	05CH-0061
Piston Guide Assembly	1	243045	05CH-0071
Plunger, for 100A/110A/110B/112/112M	1	243053	05CH-0081
Plunger Seal, UHMW-PE, for 100A/110A/110B	1	887138	05CH-0091
Plunger Seal, Black Teflon, for 100A/110A/110B	1	887138	05CH-0101
Plunger Seal, UHMW-PE, for 105/110 PREP	1	887139	05CH-0111
Plunger Seal, Black Teflon, for 105/110 PREP	1	887139	05CH-0121
Plunger Seal, UHMW-PE, for 112/112M	1	236797	05CH-0131
Plunger Seal, Black Teflon, for 112/112M	1	236797	05CH-0141
Plunger, for 114/114M/116/126	1	240714	05CH-0151
Plunger Seal, UHMW-PE, for 114/114M/116/126	1	241037	05CH-0161
Plunger Seal, Black Teflon, for 114/114M/116/126	1	241037	05CH-0171
Deuterium Lamp	1	596791	05CH-0181

- Black Teflon seals are PTFE with a carbon fill, this seal is a very good general purpose seal for any solvent.
- Gold plunger seals are made of UHMW-PE (Ultra High Molecular Weight Polyethylene) with fill that gives it a gold appearance. The seal provides good seal life and does not shed particles like that of other seals. The gold seal is inert to all common HPLC solvent mixes.