

# HPLC Instrument Parts

## Dionex HPLC Systems



05DH-0011



05DH-0087



05DH-0097



05DH-0031



05DH-0061



05DH-0101

### Ordering Information

Description	QTY	Dionex No.	Part No.
<b>Replacement parts for Dionex HPLC Systems</b>			
Inlet check valve assembly	1	38273	05DH-0011
Replacement inlet check valve cartridge	1		05DH-0021
Outlet check valve assembly	1	42761	05DH-0031
Replacement outlet check valve cartridge	1		05DH-0041
Replacement 2µm Titanium filters	10		05DH-0057
Ruby sapphire piston	1	36904	05DH-0061
UHMW-PE seals (Elastomer O-ring)	10	35686	05DH-0077
UHMW-PE seals (Hastelloy C Spring)	10	40859	05DH-0087
Backup seals	10	36901	05DH-0097
Rotor Seal (RheBuild®) kit for Dionex Ultimate 3000	1	6722.9012	05DH-0141
D2 lamp for DSA-1, VDM-1	1		05DH-0101
D2 lamp for AD20, PDA20	1		05DH-0111
D2 for UVD160, 170S, 170U, 320, 340S, 340U	1		05DH-0121
D2 lamp for Dionex VWD3000	1	6074.1110	05DH-0131

- Black Teflon seals are PTFE with a carbon fill, this seal is a very good general purpose seal for any solvent.
- Gold plunger seals are made of UHMW-PE (Ultra High Molecular Weight Polyethylene) with fill that gives it a gold appearance. The seal provides good seal life and does not shed particles like that of other seals. The gold seal is inert to all common HPLC solvent mixes.

## Hitachi/Merck HPLC Systems



05EH-0021



05EH-0061



05EH-0081



05EH-0071



05EH-0031



05EH-0112



05EH-0071



05EH-0111



05EH-0041



05EH-0051

### Ordering Information

Description	QTY	Hitachi No.	Part No.
<b>Replacement parts for Hitachi/Merck HPLC Systems</b>			
SS Check Valve Cartridges, 1/16" Ball, for L-6000/6200	2	AN0-085	05EH-0012
Inlet Check Valve Assembly, for L-6000/6200	1	AN0-0833	05EH-0021
Outlet Check Valve Assembly, for L-6000/6200	1	AN0-0834	05EH-0031
Inlet Check Valve Assembly, for L-7100	1	AN0-0836	05EH-0041
Outlet Check Valve Assembly, for L-7100	1	AN0-0837	05EH-0051
Plunger Assembly	1	810-1033	05EH-0061
Plunger Seal, Black Teflon	1	655-1080	05EH-0071
Plunger Seal, UHMW-PE	1	655-1080	05EH-0081
Plunger Seal, UHMW-PE with Hastelloy C Spring	1	655-1080	05EH-0091
Rotor Assembly for Dilutor	1	810-3085	05EH-0101
Rotor Seal Kit	1	AN0-0818	05EH-0111
Deuterium Lamp, Prealigned	1	AN0-0424	05EH-0121
Deuterium Lamp, Prealigned for Hitachi-Merck L7400/L7420/L7450/L7455/L2400	1	885-3570	05EH-0131