


HPLC Instrument Parts

Shimadzu HPLC Systems



Ordering Information			
Description	QTY	Shimadzu No.	Part No.
Replacement parts for Shimadzu HPLC Systems			
Inlet check valve assembly for LC-10AT/10ATVP	1	228-32166-91	05HH-0011
Outlet check valve assembly for LC-10AT/10ATVP	1	228-32531-91	05HH-0021
Replacement check valve cartridge, 2/Pk	2		05HH-0032
Inlet check valve assembly for LC-10ADVP/20AD/20AB	1	228-39093-92	05HH-0041
Outlet check valve assembly for LC-10ADVP/20AD/20AB	1	228-34976-91	05HH-0051
Replacement check valve cartridge, 2/Pk	2		05HH-0062
Inlet check valve assembly for LC-10AD/600/9A	1	228-18522-91	05HH-0071
Outlet check valve assembly for LC-10AD/600/9A	1	228-18522-92	05HH-0081
Replacement check valve cartridge, 2/Pk	2		05HH-0092
Inlet check valve assembly for LC-10AS/6A	1	228-12353-91	05HH-0101
Outlet check valve assembly for LC-10AS/6A	1	228-09054-93	05HH-0111
Replacement valve cartridge, 2/Pk	2		05HH-0122
Inlet check valve assembly for LC-8A	1	228-20256-91	05HH-0131
Outlet check valve assembly for LC-8A	1	228-09054-93	05HH-0141
Inlet check valve assembly for Prominence LC-20AD/B	1	228-48249-91, 228-45704-91	05HH-0371
Replacement check valve cartridge, 2/Pk	2		05HH-0152
PEEK Frit for Inline Filter for LC-10AD/ATvp/20AD/AT/AB	1	228-48607-91	05HH-0391
Sapphire Piston for LC-10ADVP/2010	1	228-35601-91	05HH-0161
Seals, UHMW-PE for LC-10ADVP/2010	1	228-32628-00	05HH-0171
Seals, Black Teflon for LC-10ADVP/2010	1	228-35146-00	05HH-0181
Sapphire Piston for LC-6A	1	228-12904-93	05HH-0191
Sapphire Piston for LC-10AS	1	228-17019-93	05HH-0201
Sapphire Piston for LC-10AT	1	228-32654-91	05HH-0211
Sapphire Piston for LC-10ATVP	1	228-35009-92	05HH-0221
Plunger assembly for LC-10ADVP/2010	1	228-35601-93	05HH-0381
Seals, UHMW-PE for LC-6A/7A/10AS/10AT/10ATVP	1	228-21975-00	05HH-0231
Seals, Black Teflon for LC-6A/7A/10AS/10AT/10ATVP	1	228-11999-00	05HH-0241
Sapphire Piston for LC-600/9A/10AD	1	228-18523-91	05HH-0251
Seals, UHMW-PE for LC-6A/7A/10AS/10AT/10ATVP	1	228-18745-00	05HH-0261
Seals, Black Teflon for LC-6A/7A/10AS/10AT/10ATVP	1	228-35146-00	05HH-0271
Needle Seal for SIL-10A/10XL/10ADVP	1	228-33355-04	05HH-0281
Rotor Seal for SIL-10A/10XL/10ADVP	1	228-21217-97	05HH-0291
Rotor Seal Assembly for SIL-10A/10XL/10ADVP	1	228-21217-91	05HH-0301
Stator Assembly for SIL-10A/10XL/10ADVP	1	228-21220-91	05HH-0311
D2 lamp for SPD-10A/10AV/M10A	1	228-34016-02	05HH-0321
D2 Lamp for SPD-2A/6A	1	062-65056-03	05HH-0331
D2 Lamp for UV-1700	1	062-65055-05	05HH-0351
Xenon Lamp for RF-10A, 10AD, 10AXL, 535, 540, 551, 1501, 5000, 5301	1	228-34216-00	05HH-0341
Tungsten lamp for SPD-10A/10AV/M10A	1	228-34110-91	05HH-0361

- Black Teflon seals are PTFE with a carbon fill, this seal is a very good general purpose seal for any solvent.
- Gold plunger seals are made of UHMW-PE (Ultra High Molecular Weight Polyethylene) with fill that gives it a gold appearance. The seal provides good seal life and does not shed particles like that of other seals. The gold seal is inert to all common HPLC solvent mixes.