

Membrane Filters and Accessories

Membrane Filters

- Variety of Pore Sizes and Membrane Types
- Non-Sterile



Hydrophilic : Aqueous and organic solvent

NYLON Membrane Filters

Size	Pore Size	QTY.	Part No.
13mm	0.2 µm	100	0231-0101
	0.45 µm	100	0231-0301
25mm	0.2 µm	100	0233-0101
	0.45 µm	100	0233-0301
47mm	0.2 µm	100	0235-0101
	0.45 µm	100	0235-0301

CA Membrane Filters

Size	Pore Size	QTY.	Part No.
13mm	0.2 µm	100	0231-1101
	0.45 µm	100	0231-1301
25mm	0.2 µm	100	0233-1101
	0.45 µm	100	0233-1301
47mm	0.2 µm	100	0235-1101
	0.45 µm	100	0235-1301

MCE Membrane Filters

Size	Pore Size	QTY.	Part No.
13mm	0.2 µm	100	0231-8101
	0.45 µm	100	0231-8301
25mm	0.2 µm	100	0233-8101
	0.45 µm	100	0233-8301
47mm	0.2 µm	100	0235-8101
	0.45 µm	100	0235-8301

Hydrophobic : Non-Aqueous, organic and high solvent resistant

PTFE Membrane Filters

Size	Pore Size	QTY.	Part No.
13mm	0.2 µm	100	0231-2101
	0.45 µm	100	0231-2301
25mm	0.2 µm	100	0233-2101
	0.45 µm	100	0233-2301
47mm	0.2 µm	100	0235-2101
	0.45 µm	100	0235-2301

PVDF Membrane Filters

Size	Pore Size	QTY.	Part No.
13mm	0.2 µm	100	0231-9101
	0.45 µm	100	0231-9301
25mm	0.2 µm	100	0233-9101
	0.45 µm	100	0233-9301
47mm	0.2 µm	100	0235-9101
	0.45 µm	100	0235-9301

Syringe Filter Holder

Use this reusable syringe filter housing with disposable filter membranes to filter samples or small quantities of mobile phases. The stainless steel or Polypropylene filter housing holds your choice of membrane. The outlet accepts any standard luer-hub needle.



Syringe Filter Holder

Description	QTY.	Part No.
13mm Stainless Steel Syringe Filter Holder	1	0210-0001
25mm Stainless Steel Syringe Filter Holder	1	0210-0002
13mm Polypropylene Syringe Filter Holder	10	0210-0003
25mm Polypropylene Syringe Filter Holder	12	0210-0004

Membrane Filters and Accessories

Glass Syringes

- Reinforced flanges with two flat sides to prevent rolling
- All borosilicate glass plunger and barrel for accurate, leak-proof performance and resistance to chemicals and thermal shock.
- Graduations are fused on to the syringe barrel
- Interchangeable plungers



5 mL Glass syringe with glass luer slip tip



20 mL Glass syringe with metal luer lock tip

Plastic Syringes

- Made of medical-grade polypropylene
- Sterile Disposable
- Ideal for use with syringe filters and general purpose applications
- Come in an easy-to-open sealed package
- Eliminate Ion Contamination in IC Applications



Glass Syringe, Luer-Tip		
Size	QTY.	Part No.
2mL	1	0211-1211
	12	0211-1213
5mL	1	0211-1411
	12	0211-1413
10mL	1	0211-1511
	12	0211-1513
20mL	1	0211-1611
	12	0211-1613
50mL	1	0211-1811
	12	0211-1813

Glass Syringe, Luer-Lock		
Size	QTY.	Part No.
2mL	1	0211-2211
	12	0211-2213
5mL	1	0211-2411
	12	0211-2413
10mL	1	0211-2511
	12	0211-2513
20mL	1	0211-2611
	12	0211-2613
50mL	1	0211-2811
	12	0211-2813

Plastic Syringe, Luer-Tip		
Size	QTY.	Part No.
1mL	100	0211-1125
3mL	100	0211-1325
5mL	100	0211-1425
10mL	100	0211-1525
20mL	50	0211-1624

Plastic Syringe, Luer-Lock		
Size	QTY.	Part No.
3mL	100	0211-2325
5mL	100	0211-2425
10mL	100	0211-2525
20mL	50	0211-2624

Filtration