

VertiPak™ Custom SPE

VertiPak™ Custom SPE can be made to order to your specification for special applications. Please specify packing type, bed weight tube volume and quantity/pack.

Ordering Instructions Example			
Requirement : Packing type : 3-Layer Carbograph/Florisil/C18			
Bed weight : 500/250/250mg			
Tube Volume : 6mL			
Quantity/PK : 30			
Packings	1 st Layer	2 nd Layer	3 rd Layer
1. Base (See page)	Carbograph	Florisil-PR	Silica
2. Particle Size (µm)	35-200	150-250	60-120
3. Pore size (Å)	-	60-80	120-140
4. Surface area (m ² /g)	100-200	300-450	300-400
5. Phase type (see page 275)	-	-	C18-LP
6. End capping (Yes or No)	-	-	Yes
7. Carbon Loading (%)	-	-	17
8. Nitrogen Loading (%)	-	-	-
9. Hydrogen Loading (%)	-	-	-
10. Bonding	-	-	Trifunctional
11. Surface coverage (µmol/m ²)	-	-	3
12. Ion-exchange capacity (meq/g)	-	-	-
13. Bed weight (mg)	500	250	250
Tubing specification			
1. Tube format (see page 276)	Tube		
2. Tube volume (see page 276)	6mL		
3. Quantity per pack	30		

Ordering Instructions			
Requirement : Packing type :			
Bed weight :			
Tube Volume :			
Quantity/PK :			
Packings	1 st Layer	2 nd Layer	3 rd Layer
1. Base (See page)			
2. Particle Size (µm)			
3. Pore size (Å)			
4. Surface area (m ² /g)			
5. Phase type (see page 275)			
6. End capping (Yes or No)			
7. Carbon Loading (%)			
8. Nitrogen Loading (%)			
9. Hydrogen Loading (%)			
10. Bonding			
11. Surface coverage (µmol/m ²)			
12. Ion-exchange capacity (meq/g)			
13. Bed weight (mg)			
Tubing specification			
1. Tube format (see page 276)			
2. Tube volume (see page 276)			
3. Quantity per pack			