


12x32mm Snap/Crimp Wide Mouth Vials

12x32mm Snap/Crimp Wide Mouth Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass or high quality polypropylene
- Total volume of 2mL and usable volume of 1.5mL
- Glass vials are available with or without spot
- Wide mouth with Snap/Crimp neck
- Aluminum or Polypropylene open-top crimp seals with choices of septa
- Limited volume inserts are available for small volume samples
- Pre-assembled vial kits include vials, caps with septa and a PP gridded box
- Compatible with Agilent, Waters®, HTA, Shimadzu, Thermo and Varian autosampler



Ordering Information

Description	QTY	Part No.
12x32mm Snap/Crimp Wide Mouth Vials		
Clear glass	100	0705-1011
Clear glass, spot	100	0705-1021
Amber glass	100	0705-1041
Amber glass, spot	100	0705-1051
Crimp Seals with Septa		
Al seals, PTFE/Red Rubber	100	0703-4011
Al seals, PTFE/White Silicone	100	0703-4021
Al seals, PTFE/Silicone/PTFE	100	0703-4031
Al seals, PTFE/Silicone, Slit	100	0703-4041
PP clear seals, PTFE/Red Rubber	100	0705-4011
PP clear seals, PTFE/White Silicone	100	0705-4021
PP clear seals, PTFE/Silicone/PTFE	100	0705-4031
PP clear seals, PTFE/Silicone, Slit	100	0705-4041
Vial Kits in Gridded Box		
Clear glass vials with PP seals and PTFE/Red Rubber Septa	100	0705-5011
Amber glass vials with PP seals and PTFE/Red Rubber Septa	100	0705-5021
Clear glass vials with PP seals and PTFE/White Silicone Septa	100	0705-5031
Amber glass vials with PP seals and PTFE/White Silicone Septa	100	0705-5041
Inserts		
Glass, 100µL, w/Bottom Spring	100	0701-6011
Glass, 250µL, w/Bottom Spring	100	0701-6021
Glass, 300µL	100	0701-6031
Glass, 400µL, Flat Bottom	100	0701-6041