


# 12x32mm Crimp Top Standard Mouth Vials

## 12x32mm Crimp Top Standard Mouth Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass or high quality polypropylene
- Total volume of 2mL and useable volume 1.5mL
- Glass vials are available with or without spot
- Standard mouth with Crimp Top neck
- Aluminum or Polypropylene open-top crimp caps with choices of septa
- Limited volume inserts are available for small volume samples
- A vial kit in PP Gridded Box includes vials, Caps with septa
- Compatible with Agilent, Waters®, HTA, Shimadzu, Thermo and Varian autosampler



### Ordering Information

Description	QTY	Part No.
<b>12x32mm Crimp Top Standard Mouth Vials</b>		
Clear glass, spot	100	0703-1011
Amber glass, spot	100	0703-1051
<b>Crimp Seals with Septa</b>		
Al seals, PTFE/Red Rubber	100	0703-4011
Al seals, PTFE/White Silicone	100	0703-4021
Al seals, PTFE/Silicone/PTFE	100	0703-4031
Al seals, PTFE/Silicone, Slit	100	0703-4041
<b>Vial Kits in Gridded Box</b>		
Clear glass spot vials with Al seals and PTFE/Red Rubber Septa	100	0703-5011
Amber glass spot vials with Al seals and PTFE/Red Rubber Septa	100	0703-5021
Clear glass spot vials with Al seals and PTFE/White Silicone Septa	100	0703-5031
Amber glass spot vials with Al seals and PTFE/White Silicone Septa	100	0703-5041
<b>Inserts</b>		
Glass, 50µL, w/Bottom Spring	100	0700-6011
Glass, 100µL, w/Bottom Spring	100	0700-6021
Glass, 150µL	100	0700-6031
Glass, 250µL, Flat Bottom	100	0700-6041