


# 12x32mm Screw Thread Wide Mouth Vials

## 12x32mm Screw Thread Wide Mouth Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass or high quality polypropylene
- Total volume of 2mL and usable volume of 1.5mL
- Glass vials are available with or without spot
- Wide mouth with 10/425 Screw Thread
- Polypropylene open-top screw caps with choices of septa
- Limited volume inserts are available for small volume samples
- A vial kit in PP Gridded Box includes vials, Caps with septa
- Compatible with Waters®, Shimadzu, Perkin Elmer, Jasco and Varian autosampler



### Ordering Information

Description	QTY	Part No.
<b>12x32mm Screw Thread Wide Mouth Vials</b>		
Clear glass	100	0702-1011
Clear glass, spot	100	0702-1021
Clear glass, silanized, spot	100	0702-1031
Amber glass	100	0702-1041
Amber glass, spot	100	0702-1051
Amber glass, silanized, spot	100	0702-1061
PP	100	0702-1071
<b>10/425 Open-Top PP Screw Caps</b>		
Black	100	0702-2011
Red	100	0702-2021
<b>Septa for 10/425 Caps</b>		
PTFE/Red Rubber, 0.04" Thickness	100	0702-3011
PTFE/White Silicone, 0.06" Thickness	100	0702-3021
PTFE/Silicone/PTFE, 0.06" Thickness	100	0702-3031
PTFE/Silicone, 0.06" Thickness, with Slit	100	0702-3041
PTFE, 0.01" Thickness	100	0702-3051
<b>Caps with Septa</b>		
Black open-top caps, PTFE/Red Rubber	100	0702-4011
Black open-top caps, PTFE/White Silicone	100	0702-4021
Black open-top caps, PTFE/Silicone/PTFE	100	0702-4031
Black open-top caps, PTFE/Silicone, Slit	100	0702-4041
White Solid-top phenolic caps, PTFE/silicone	100	0702-4051
Clear Solid-top integral caps with septa	100	0702-4061
<b>Vial Kits in Gridded Box</b>		
Clear glass spot vials with black caps and PTFE/White Silicone Septa	100	0702-5011
Amber glass spot vials with black caps and PTFE/White Silicone Septa	100	0702-5021
<b>Inserts</b>		
Glass, 100µL, w/Bottom Spring	100	0701-6011
Glass, 250 µL, w/Bottom Spring	100	0701-6021
Glass, 300 µL	100	0701-6031
Glass, 400 µL, Flat Bottom	100	0701-6041