


# 12x32mm Screw Thread Standard Mouth Vials

## 12x32mm Screw Thread Standard Mouth Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass or high quality polypropylene
- Total volume of 2 mL and usable volume of 1.5mL
- Glass vials are available with or without spot
- Standard mouth with 8/425 Screw Thread
- Polypropylene open-top screw caps with choices of septa
- Limited volume inserts are available for small volume samples
- A vial kit in PP Gridded Box includes vials, Caps with septa
- Compatible with Waters® Alliance®, Shimadzu, Hitachi, Jasco and Varian autosampler



### Ordering Information

Description	QTY	Part No.
<b>12x32mm Screw Thread Standard Mouth Vials</b>		
Clear glass	100	0700-1011
Clear glass, spot	100	0700-1021
Clear glass, silanized, spot	100	0700-1031
Amber glass	100	0700-1041
Amber glass, spot	100	0700-1051
Amber glass, silanized, spot	100	0700-1061
PP	100	0700-1071
<b>8/425 Open-Top PP Caps</b>		
Black	100	0700-2011
Red	100	0700-2021
Green	100	0700-2031
Yellow	100	0700-2041
<b>Septa for 8/425 Caps</b>		
PTFE/Red Rubber, 0.04" Thickness	100	0700-3011
PTFE/White Silicone, 0.06" Thickness	100	0700-3021
PTFE/Silicone/PTFE, 0.06" Thickness	100	0700-3031
PTFE/Silicone, 0.06" Thickness with Slit	100	0700-3041
PTFE, 0.01" Thickness	100	0700-3051
<b>Caps with Septa</b>		
Black Open-Top Cap, PTFE/Red Rubber	100	0700-4011
Black Open-Top Cap, PTFE/White Silicone	100	0700-4021
Black Open-Top Cap, PTFE/Silicone/PTFE	100	0700-4031
Black Open-Top Cap, PTFE/Silicone, Slit	100	0700-4041
Yellow Open-Top Cap, PTFE/Silicone	100	0700-4051
Black Solid-Top Cap, PTFE/Silicone	100	0700-4061
<b>Vial Kits in Gridded Box</b>		
Clear glass spot vials with black open-top caps and PTFE/White Silicone Septa	100	0700-5011
Amber glass spot vials with black open-top caps and PTFE/White Silicone Septa	100	0700-5021
Clear glass spot vials with yellow open-top caps and PTFE/White Silicone Septa	100	0700-5031
Amber glass spot vials with yellow open-top caps and PTFE/White Silicone Septa	100	0700-5041
<b>Inserts</b>		
Glass, 50µL, w/Bottom Spring	100	0700-6011
Glass, 100µL, w/Bottom Spring	100	0700-6021
Glass, 150µL	100	0700-6031
Glass, 250µL, Flat Bottom	100	0700-6041