


# 15x45mm Screw Thread Vials

## 15x45mm Screw Thread Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass or high quality polypropylene
- Total volume of 5.2mL and useable volume of 4mL
- Glass vials are available with or without spot
- Wide mouth with 13/425 Screw Thread
- Polypropylene open-top screw caps with choices of septa
- Limited volume inserts are available for small volume samples
- Pre-assembled vial kits include vials and Caps with septa
- Compatible with Waters®, Shimadzu, Perkin Elmer, Jasco and Varian autosampler



### Ordering Information

Description	QTY	Part No.
<b>15x45mm Screw Thread Vials</b>		
Clear glass	100	0706-1011
Clear glass, spot	100	0706-1021
Clear glass, silanized, spot	100	0706-1031
Amber glass	100	0706-1041
Amber glass, spot	100	0706-1051
Amber glass, silanized, spot	100	0706-1061
PP, Well Base, 2mL	100	0706-1071
<b>13/425 Screw Caps</b>		
Black, PP, Open-Top	100	0706-2011
White, PP, Open-Top	100	0706-2021
<b>Septa for 13/425 Screw Caps</b>		
PTFE/Red Rubber	100	0706-3011
PTFE/White Silicone	100	0706-3021
PTFE/Silicone/PTFE	100	0706-3031
PTFE/Silicone with Slit	100	0706-3041
PTFE 10mil	100	0706-3041
<b>Caps with Septa</b>		
Black Open-Top PP Caps, PTFE/Red Rubber	100	0706-4011
Black Open-Top PP Caps, PTFE/Silicone	100	0706-4021
Black Open-Top PP Caps, PTFE/Silicone/PTFE	100	0706-4031
Black Open-Top Caps, PTFE/Silicone	100	0706-4041
Black Solid-Top Phenolic Caps, PTFE/Silicone	100	0706-4051
<b>Vial Kits</b>		
Clear glass spot vials with black caps and PTFE/White Silicone Septa	100	0706-5011
Amber glass spot vials with black caps and PTFE/White Silicone Septa	100	0706-5021
<b>Inserts</b>		
Glass, 300µL, Flanged	100	0706-6011
PP, 300µL, Flanged	100	0706-6021
Glass, 750µL, w/Bottom Spring	100	0706-6031
Glass, 1mL, Flat Bottom	100	0706-6041