


15x45mm Snap/Crimp Vials

15x45mm Snap/Crimp Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass or high quality polypropylene
- Total volume of 5.2mL and useable volume of 4mL
- Glass vials are available with or without spot
- Wide mouth with snap/crimp neck style
- Snap/crimp seals with choices of septa
- Limited volume inserts are available for small volume samples
- Pre-assembled vial kits include vials and Caps with septa
- Compatible with Waters® WISP™ 48 Style and Shimadzu autosampler



Ordering Information

Description	QTY	Part No.
15x45mm Snap/Crimp Vials		
Clear glass	100	0707-1011
Clear glass, spot	100	0707-1021
Amber glass	100	0707-1041
Amber glass, spot	100	0707-1051
PP, Well Base, 2mL	100	0707-1071
Snap/Crimp Seals with Septa		
Al seals, PTFE/Red Rubber	100	0707-4011
Al seals, PTFE/White Silicone	100	0707-4021
Al seals, PTFE/Silicone/PTFE	100	0707-4031
Al seals, PTFE/Silicone, Slit	100	0707-4041
PP clear seals, PTFE/Red Rubber	100	0707-4051
PP clear seals, PTFE/White Silicone	100	0707-4061
PP clear seals, PTFE/Silicone/PTFE	100	0707-4071
PP clear seals, PTFE/Silicone, Slit	100	0707-4081
Vial Kits in Gridded Box		
Clear glass vials with Al seals and PTFE/White Silicone Septa	100	0707-5011
Amber glass vials with Al seals and PTFE/White Silicone Septa	100	0707-5021
Clear glass vials with PP seals and PTFE/White Silicone Septa	100	0707-5031
Amber glass vials with PP seals and PTFE/White Silicone Septa	100	0707-5041
Inserts		
Glass, 300µL, Flanged	100	0706-6011
PP, 300µL, Flanged	100	0706-6021
Glass, 750µL, w/Bottom Spring	100	0706-6031
Glass, 1mL, Flat Bottom	100	0706-6041