


# 15x45mm Screw Thread Vials

## 15x45mm Screw Thread Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass or high quality polypropylene
- Total volume of 5.2mL and useable volume of 4mL
- Glass vials are available with or without spot
- Wide mouth with 13/425 Screw Thread
- Polypropylene open-top screw caps with choices of septa
- Limited volume inserts are available for small volume samples
- Pre-assembled vial kits include vials and Caps with septa
- Compatible with Waters®, Shimadzu, Perkin Elmer, Jasco and Varian autosampler



### Ordering Information

| Description                                                          | QTY | Part No.  |
|----------------------------------------------------------------------|-----|-----------|
| <b>15x45mm Screw Thread Vials</b>                                    |     |           |
| Clear glass                                                          | 100 | 0706-1011 |
| Clear glass, spot                                                    | 100 | 0706-1021 |
| Clear glass, silanized, spot                                         | 100 | 0706-1031 |
| Amber glass                                                          | 100 | 0706-1041 |
| Amber glass, spot                                                    | 100 | 0706-1051 |
| Amber glass, silanized, spot                                         | 100 | 0706-1061 |
| PP, Well Base, 2mL                                                   | 100 | 0706-1071 |
| <b>13/425 Screw Caps</b>                                             |     |           |
| Black, PP, Open-Top                                                  | 100 | 0706-2011 |
| White, PP, Open-Top                                                  | 100 | 0706-2021 |
| <b>Septa for 13/425 Screw Caps</b>                                   |     |           |
| PTFE/Red Rubber                                                      | 100 | 0706-3011 |
| PTFE/White Silicone                                                  | 100 | 0706-3021 |
| PTFE/Silicone/PTFE                                                   | 100 | 0706-3031 |
| PTFE/Silicone with Slit                                              | 100 | 0706-3041 |
| PTFE 10mil                                                           | 100 | 0706-3041 |
| <b>Caps with Septa</b>                                               |     |           |
| Black Open-Top PP Caps, PTFE/Red Rubber                              | 100 | 0706-4011 |
| Black Open-Top PP Caps, PTFE/Silicone                                | 100 | 0706-4021 |
| Black Open-Top PP Caps, PTFE/Silicone/PTFE                           | 100 | 0706-4031 |
| Black Open-Top Caps, PTFE/Silicone                                   | 100 | 0706-4041 |
| Black Solid-Top Phenolic Caps, PTFE/Silicone                         | 100 | 0706-4051 |
| <b>Vial Kits</b>                                                     |     |           |
| Clear glass spot vials with black caps and PTFE/White Silicone Septa | 100 | 0706-5011 |
| Amber glass spot vials with black caps and PTFE/White Silicone Septa | 100 | 0706-5021 |
| <b>Inserts</b>                                                       |     |           |
| Glass, 300µL, Flanged                                                | 100 | 0706-6011 |
| PP, 300µL, Flanged                                                   | 100 | 0706-6021 |
| Glass, 750µL, w/Bottom Spring                                        | 100 | 0706-6031 |
| Glass, 1mL, Flat Bottom                                              | 100 | 0706-6041 |