

Crimpers and Decappers

Electronic Crimper and Decapper

- Weighs less than 1-1/2lbs
- Crimps hundreds of crimp top vials per battery charge
- Easy to use, with a comfortable, ergonomic grip
- Consistent and reproducible crimping
- Settings are saved when not in use
- Available for use with crimp styles vials sizes of 11mm, 13mm and 20mm
- Provided with universal power supply/recharger



0716-7020

Plier Decapper

- Solid stainless construction
- Low cost hand tool which allows crimped aluminum seals to be removed quickly and easily from vials and bottles.
- Available for use with crimp styles vials sizes of 8mm, 11mm, 13mm and 20mm



0716-7170

Hand Crimper and Decapper

- Solid construction, combined with quality components manufactured to precise tolerances, provide consistent, dependable performance.
- Built-in adjustable stop allows the user to repeatedly achieve an optimum seal.
- Use to crimp or decap vials and bottles in one quick, easy motion
- Available for use with crimp styles vials sizes of 8mm, 11mm, 13mm and 20mm



0716-7080



0716-7120

Ordering Information		
Description	QTY	Part No.
Electronic Crimpers for Crimp Top Vials		
11mm Electronic Crimper	1	0716-7010
13mm Electronic Crimper	1	0716-7020
20mm Electronic Crimper	1	0716-7030
Electronic Decappers for Crimp Top Vials		
11mm Electronic Decapper	1	0716-7040
13mm Electronic Decapper	1	0716-7050
20mm Electronic Decapper	1	0716-7060
Hand Crimpers for Crimp Top Vials		
8mm Hand Crimper	1	0716-7070
11mm Hand Crimper	1	0716-7080
13mm Hand Crimper	1	0716-7090
20mm Hand Crimper	1	0716-7100
Hand Decappers for Crimp Top Vials		
8mm Hand Decapper	1	0716-7110
11mm Hand Decapper	1	0716-7120
13mm Hand Decapper	1	0716-7130
20mm Hand Decapper	1	0716-7140
Plier Decappers for Crimp Top Vials		
8mm Plier Decapper	1	0716-7150
11mm Plier Decapper	1	0716-7160
13mm Plier Decapper	1	0716-7170
20mm Plier Decapper	1	0716-7180