

Storage Vials



- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass with screw thread finish
- Choice of 2mL, 4mL, 8mL, 12mL, 16mL and 24mL
- Solid Top caps are phenolic or PP with liners
- Caps are autoclavable

Ordering Information

| Description | QTY | Part No. |
|--|-----|-----------|
| Screw Thread Storage Vials | | |
| 2mL Clear Glass, 12x32mm, 8/425 Finish | 100 | 0700-1011 |
| 2mL Amber Glass, 12x32 mm, 8/425 Finish | 100 | 0700-1041 |
| 4mL Clear Glass, 15x45mm, 13/425 Finish | 100 | 0706-1011 |
| 4mL Amber Glass, 15x45mm, 13/425 Finish | 100 | 0706-1041 |
| 8mL Clear Glass, 17x60mm, 15/425 Finish | 100 | 0711-1011 |
| 8mL Amber Glass, 17x60mm, 15/425 Finish | 100 | 0711-1021 |
| 12mL Clear Glass, 19x65mm, 15/425 Finish | 100 | 0711-1031 |
| 12mL Amber Glass, 19x65mm, 15/425 Finish | 100 | 0711-1041 |
| 16mL Clear Glass, 21x70mm, 18/400 Finish | 100 | 0711-1051 |
| 16mL Amber Glass, 21x70mm, 18/400 Finish | 100 | 0711-1061 |
| 24mL Clear Glass, 23x85mm, 20/400 Finish | 100 | 0711-1071 |
| 24mL Amber Glass, 23x85mm, 20/400 Finish | 100 | 0711-1081 |
| Solid Top Screw Caps with Liners | | |
| PP, 8/425 with PTFE Liners | 100 | 0711-4011 |
| PP, 13/425 with PTFE Liners | 100 | 0711-4021 |
| PP, 15/425 with PTFE Liners | 100 | 0711-4031 |
| PP, 18/400 with PTFE Liners | 100 | 0711-4041 |
| PP, 20/400 with PTFE Liners | 100 | 0711-4051 |
| Phenolic, 8/425 with PTFE/F217 Liners | 100 | 0711-4061 |
| Phenolic, 13/425 with PTFE/F217 Liners | 100 | 0711-4071 |
| Phenolic, 15/425 with PTFE/F217 Liners | 100 | 0711-4081 |
| Phenolic, 18/400 with PTFE/F217 Liners | 100 | 0711-4091 |
| Phenolic, 20/400 with PTFE/F217 Liners | 100 | 0711-4101 |

Culture Tubes



- Round-Bottom Screw-Top Tubes, Type I Class A Borosilicate glass
- Phenolic Cap with Teflon-Faced Rubber or White Rubber
- For tissue culture and general bacteriology laboratory use

Ordering Information

| Description | QTY | Part No. |
|---|-----|-----------|
| Round-Bottom Screw-Top Culture Tubes | | |
| 8mL, Clear Glass, 13x100mm, 13/415 Finish | 250 | 0714-1012 |
| 10mL, Clear Glass, 16x100mm, 15/415 Finish | 250 | 0714-1022 |
| 15mL, Clear Glass, 16x125mm, 15/415 Finish | 250 | 0714-1032 |
| Phenolic Cap | | |
| Cap with Teflon-Faced Rubber, 13/415 Finish | 250 | 0714-2012 |
| Cap with Teflon-Faced Rubber, 15/415 Finish | 250 | 0714-2022 |
| Cap with White Rubber, 15/415 Finish | 250 | 0714-2032 |