


Crimp Top Headspace Vials

Crimp Top Headspace Vials

- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass with standard or tapered 20mm finish
- Choice of 6mL, 10 mL and 20mL
- Round bottom vials are compatible with most autosampler
- Flat bottom vials are preferred for manual headspace applications
- Pre-assembled vial kits include vials and seals with septa
- Compatible with Waters®, Shimadzu, PerkinElmer, Jasco and Varian autosampler



Ordering Information

Description	QTY	Part No.
Crimp Top Headspace Vials		
6mL Clear Glass, Flat Bottom, Standard	100	0708-1011
6mL Clear Glass, Flat Bottom, Tapered	100	0708-1021
6mL Clear Glass, Round Bottom, Tapered	100	0708-1031
10mL Clear Glass, Flat Bottom, Standard	100	0708-1041
10mL Clear Glass, Flat Bottom, Tapered	100	0708-1051
10mL Clear Glass, Round Bottom, Tapered	100	0708-1061
10mL Amber Glass, Round Bottom, Tapered	100	0708-1071
20mL Clear Glass, Flat Bottom, Standard	100	0708-1081
20mL Clear Glass, Flat Bottom, Tapered	100	0708-1091
20mL Clear Glass, Round Bottom, Tapered	100	0708-1101
20mL Amber Glass, Round Bottom, Tapered	100	0708-1111
20mm Crimp Seals		
Al, Center Hole	100	0708-2011
Al, Center Pressure Release	100	0708-2021
Al, Center Tear-Off	100	0708-2031
Al, Center Tear-Away	100	0708-2041
Al, Center Hole, Magnetic	100	0708-2051
20mm Septa and Stopper		
PTFE/Gray Butyl	100	0708-3011
PTFE/Red Rubber	100	0708-3021
PTFE/White Silicone	100	0708-3031
PTFE/Clear Silicone	100	0708-3041
PTFE/Gray Butyl, Molded Pharmafix	100	0708-3051
Red Rubber Stopper	100	0708-3061
Gray Butyl Stopper	100	0708-3071
20mm Crimp Seals with Septa		
Al, Center Hole with PTFE/Red Rubber	100	0708-4011
Al, Center Hole with PTFE/White Silicone	100	0708-4021
Al, Center Pressure Release with PTFE/White Silicone	100	0708-4031
Al, Gold Magnetic with PTFE/White Silicone	100	0708-4041