

# Vial Racks

## Stackable Racks

- Ideal for storing and transporting 8-28mm diameter chromatography vials
- 50-96 indexed positions for easy sample identification
- Stackable and autoclavable polypropylene



0713-1031

## Snap Racks

- Ideal for storing and transporting 8-28mm diameter chromatography vials
- 50-96 indexed positions for easy sample identification
- Stackable and autoclavable polypropylene



0713-2031

## Storage Boxes

- Hold 100 12x32mm vials
- Clear polypropylene
- Available in Blue, Red, Yellow, Violet, Green and Natural
- Raised edge provides non-slip, stackable storage with writing panel on lid front for ID
- Hinged lid with secure snap clasp
- Temp range -70 - 140°C, ideal for freezer or room storage
- Autoclavable
- Dimension of 5.5"L x 6.0"W x 2.1"H



0713-3075

Ordering Information		
Description	QTY	Part No.
<b>Stackable Racks</b>		
8mm, 96holds, 185x127x23mm	1	0713-1011
8mm, 96holds, 185x127x23mm	5	0713-1025
12mm, 50holds, 194x100x23mm	1	0713-1031
12mm, 50holds, 194x100x23mm	5	0713-1045
15mm, 48holds, 270x97x30mm	1	0713-1051
15mm, 48holds, 270x97x30mm	5	0713-1065
15mm, 90holds, 327x172x30mm	1	0713-1071
15mm, 90holds, 327x172x30mm	5	0713-1085
28mm, 50holds, 325x171x30mm	1	0713-1091
28mm, 50holds, 325x171x30mm	5	0713-1105
<b>Snap Racks</b>		
8mm, 96holds, White	1	0713-2011
8mm, 96holds, Blue	1	0713-2021
12mm, 50holds, White	1	0713-2031
12mm, 50holds, Blue	1	0713-2041
<b>Storage Boxes</b>		
12mm, 100holds, Blue	5	0713-3015
12mm, 100holds, Red	5	0713-3025
12mm, 100holds, Violet	5	0713-3035
12mm, 100holds, Yellow	5	0713-3045
12mm, 100holds, Natural	5	0713-3055
12mm, 100holds, Green	5	0713-3065
12mm, 100holds, Assorted, 1 of EA/5 Colors	5	0713-3075