

Storage Vials



- Available in clear, Type I Class A or amber, Type I Class B Borosilicate glass with screw thread finish
- Choice of 2mL, 4mL, 8mL, 12mL, 16mL and 24mL
- Solid Top caps are phenolic or PP with liners
- Caps are autoclavable

Ordering Information

Description	QTY	Part No.
Screw Thread Storage Vials		
2mL Clear Glass, 12x32mm, 8/425 Finish	100	0700-1011
2mL Amber Glass, 12x32 mm, 8/425 Finish	100	0700-1041
4mL Clear Glass, 15x45mm, 13/425 Finish	100	0706-1011
4mL Amber Glass, 15x45mm, 13/425 Finish	100	0706-1041
8mL Clear Glass, 17x60mm, 15/425 Finish	100	0711-1011
8mL Amber Glass, 17x60mm, 15/425 Finish	100	0711-1021
12mL Clear Glass, 19x65mm, 15/425 Finish	100	0711-1031
12mL Amber Glass, 19x65mm, 15/425 Finish	100	0711-1041
16mL Clear Glass, 21x70mm, 18/400 Finish	100	0711-1051
16mL Amber Glass, 21x70mm, 18/400 Finish	100	0711-1061
24mL Clear Glass, 23x85mm, 20/400 Finish	100	0711-1071
24mL Amber Glass, 23x85mm, 20/400 Finish	100	0711-1081
Solid Top Screw Caps with Liners		
PP, 8/425 with PTFE Liners	100	0711-4011
PP, 13/425 with PTFE Liners	100	0711-4021
PP, 15/425 with PTFE Liners	100	0711-4031
PP, 18/400 with PTFE Liners	100	0711-4041
PP, 20/400 with PTFE Liners	100	0711-4051
Phenolic, 8/425 with PTFE/F217 Liners	100	0711-4061
Phenolic, 13/425 with PTFE/F217 Liners	100	0711-4071
Phenolic, 15/425 with PTFE/F217 Liners	100	0711-4081
Phenolic, 18/400 with PTFE/F217 Liners	100	0711-4091
Phenolic, 20/400 with PTFE/F217 Liners	100	0711-4101

Culture Tubes



- Round-Bottom Screw-Top Tubes, Type I Class A Borosilicate glass
- Phenolic Cap with Teflon-Faced Rubber or White Rubber
- For tissue culture and general bacteriology laboratory use

Ordering Information

Description	QTY	Part No.
Round-Bottom Screw-Top Culture Tubes		
8mL, Clear Glass, 13x100mm, 13/415 Finish	250	0714-1012
10mL, Clear Glass, 16x100mm, 15/415 Finish	250	0714-1022
15mL, Clear Glass, 16x125mm, 15/415 Finish	250	0714-1032
Phenolic Cap		
Cap with Teflon-Faced Rubber, 13/415 Finish	250	0714-2012
Cap with Teflon-Faced Rubber, 15/415 Finish	250	0714-2022
Cap with White Rubber, 15/415 Finish	250	0714-2032